

## Asepti Active Wipes

Hospital Grade Disinfectant Wipe - 325 x Pail, 100 x Soft Pack

ARTG 340994







## Features

- Asepti Active Wipes are a hospital grade disinfectant approved for single step cleaning and disinfection of dirty surfaces
- Asepti Active Wipes are a broad spectrum bactericide effective against VRE & MRSA and has virucidal efficacy against Human Coronavirus, COVID-19 [SARS-CoV-2] and Influenza A
- Asepti Active Wipes contains a combination of quaternary biocides - Benzalkonium Chloride & Dialkyldimethylammonium chloride - designed to provide optimal efficacy against QAC resistant organisms
- Asepti Active Wipes are easy to use, they ensure the correct concentration of chemistry is available in each application and reduce the risk of cross contamination
- Asepti Active Wipes are conveniently packed, low lint, low foam and are available in a range of pack sizes to suit varying applications
- Asepti Active Wipes are compatible with Asepti Active Liquid
- A range of dispensers to suit the [flat pack](#) and [pail](#) format are also available encouraging accessibility and promoting usage

### Conditions & Contact Times

- Broad spectrum bactericidal efficacy including Pseudomonas aeruginosa; Proteus vulgaris, Escherichia coli; Staphylococcus aureus: Contact Time: 8 minutes in dirty conditions
- Bacterial outbreak event, Antibiotic resistant bacteria: Vancomycin Resistant; Enterococcus faecalis - VRE; Methicillin Resistant; Staphylococcus aureus - MRSA: Contact Time: 5 minutes in dirty conditions
- Viral outbreak event: Human Coronavirus COVID-19 [SARS-CoV-2] severe acute respiratory syndrome coronavirus 2; Influenza A: Contact Time: 5 minutes in clean conditions; 10 minutes in dirty conditions

Product Code	Product Description	Pack Size	Quote
--------------	---------------------	-----------	-------

Product Code	Product Description	Pack Size	Quote
7101608	Asepti Active Wipes, Soft Pack	12 x 100 wipes	
7400573	Asepti Active Wipes, Pail	3 x 325 wipes	
99001137, 10111807	Plastic Holder to suit Flat Pack	ea	
9996-1208, 10112033	Wall Bracket to suit Pail	ea	
99960090, 10111973	Wire Holder to suit Flat Pack	ea	
HC-00023, 10112750	Procedure & Technique Wall Chart	ea	

## Safety Information

Do not mix with any other products. For further information and safety instructions see product label and Safety Data Sheet.

ALWAYS READ THE LABEL. FOLLOW DIRECTIONS OF USE.

## Directions For Use

- Wear gloves and additional PPE where required.
- Asepti Active Wipes are safe for use on most surfaces. They are not suitable for therapeutic devices or as a baby wipe.
- Asepti Active Wipes are single use only, supporting infection control protocols.
- Do not dispose of wipes in toilets.

## Related Products

- [Asepti Active Liquid](#)
- [Asepti Neutral Wipes](#)
- [Asepti Alcohol Wipes](#)

## Applications

- Patient Room
- Ward Areas
- Patient / Public Bathrooms
- Theatre
- ICU
- Emergency Department
- Central Sterile Department

## Resources

---

- SAFETY DATA SHEET (SDS)
- Video
- Brochure