





## Erlab Halo Sense

Laboratory Air Quality Sensor for Volatile Organic Compounds (VOCs), Acid and Formaldehyde, to keep your laboratory air quality safe



## Features

- Erlab Halo Sense has been designed to continuously monitor the lab air quality for a large spectrum of chemical vapors, to keep you and your laboratory personnel safe.
- Halo Sense allows you to be informed about the air you breath in the lab by offering 3 sensors to choose from; Volatile Organic Compounds (VOCs), Acid and Formaldehyde.
- When Halo Sense detects chemical vapours in the zone, a soft band of LED light blinks and beeps are emitted to warn you that the air is polluted with solvents, acids or formaldehyde.
- Erlab Halo Sense can be connected to either a computer or the Erlab eGuard app for control and monitoring.
- Fixed to a wall or placed on a table, Halo Sense helps you to assure your laboratory personnel that they are working in a safe environment.

Product Code	Product Description	Pack Size	Quote
FCC-HS-VOC	Erlab Halo Sensor, VOC Air Quality Sensor	ea	
FCC-HS-FOR	Erlab Halo Sense, Formaldehyde Air Quality Sensor	ea	
FCC-HS-ACID	Erlab Halo Sense, Acid Air Quality Sensor	ea	

## Technical Data

Dimensions (mm)	
<b>Halo Sense, VOC</b>	
Product Outer (H x W x D):	181 x 181 x 67 mm
Sensors	Volatile organic vapours
Power Consumption W	10
<b>Halo Sense, VOC</b>	

Dimensions (mm)	
Product Outer (H x W x D):	181 x 181 x 67 mm
Gas / Vapours	Formaldehyde vapours
Power Consumption W	10
<b>Halo Sense, Acid</b>	
Product Outer (H x W x D):	181 x 181 x 67 mm
Gas / Vapours	Acid vapours
Power Consumption W	10

## Relevant Product

- [Halo Air Purifier](#)
- [ChemDAQ SteriTrac 2+](#)
- [Captair Smart 321S](#)
- [Captair Smart 321](#)
- [Captair Smart 391](#)
- [Captair Smart 392](#)
- [Captair Smart 481](#)
- [Captair Smart 483](#)
- [Captair Smart 632](#)
- [Captair Smart 633](#)
- [Captair Smart 714](#)
- [Captair Bio](#)
- [Captair Flow](#)
- [Captair Pyramid](#)
- [Captair Weighing](#)
- [Green Fume Hood](#)
- [Captair Store 834](#)
- [Captair Ministore](#)

## Applications

- Life Science
- Healthcare Facilities
- Hospitals

## Resources

- Video
- Brochure