

Meiko TopClean 60

TopClean 60 is an innovative multiwasher and disinfection machine for just about anything



Features

- TopClean 60 offerings outstanding hygienic results through its dedicated Chemo-thermal disinfection, which is ideal for disinfecting heat-sensitive theatre boots
- Provides the best method of disinfection for most applications with an A0 value of 60
- TopClean 60s rapid disinfection process allows for washed items to circulate back into use quickly and efficiently
- The simple to use keypad and display provides user-friendly operation
- Various programs to suit different cleaning needs; Short cycle, Standard cycle or Intensive cycle
- TopClean 60 uses minimal precious resources such as water, electricity and chemicals, thus saving you time and money
- TopClean 60 offers the ability to change the water for tackling heavily soiled items using the intensive cycle
- Automatic door lock and unlock system, providing additional security
- The innovative technology of AktivPlus fine filters provide outstanding cleaning results
- TopClean 60s large capacity can cater for the disinfection of transport containers and 10 litre buckets
- Universal basket & trays provide optimum use of washing chamber and versatility

Product Code	Product Description	Pack Size	Quote
MWC-10-000007 34	TopClean 60 Multipurpose Washer + SK500 basket	ea	

Relevant Product

- [Asepti Synergy WD, 5 L](#)
- [Asepti Autozyme Pr, 5 L](#)
- [Asepti Synergy 5, 5L](#)
- [Asepti Wash, 5L](#)
- [Asepti Sonic, 5L](#)
- [Asepti RHW, 5L](#)
- [Asepti Rinse Clean, 5L](#)

Applications

- Infection Control
- Healthcare Facilities
- Residential Aged Care

Resources

- Video
- Video
- Brochure