


Meiko Topic 20.3

Meiko Topic machines let you focus on caring for your residents by taking care of the time consuming cleaning chores



Features

- Topic 20 is a compact hygiene solution offered in a floor standing format for easy installation with minimal space required
- Budget friendly compact solution for bed pan cleaning that packs a punch in hygiene outcomes
- Topic 20 is ideal for refurbishment and reno projects as it's quick and easy to install
- Ideal for nursing homes and aged care facilities
- Topic 20 versatility also allows you to clean other items aside from bed pans like washbowls and children's toys
- Topic 20 A model offers the inclusion of a slop sink for hygiene disposal
- Flexible integration to replace out-of-date old technology; Topic can be easily modified to use the same connections
- Offering safe and hygienic results due to value A0 control system to DIN EN IS 1588e
- Individually adjustable programme cycles for short, economical turnaround
- Topic 20 offers outstanding and thorough cleaning results with their power nozzle system and patented steam-proof cleaning & disinfection
- Thermal disinfection provides hygiene assurance along with multifunctional utensil holders for maximum flexibility
- Self cleaning wash chamber design with no corners or tight angles
- Rest assured bedpans, urine bottles and commode buckets are emptied, cleaned and hygienically disinfected, reliably in a sealed environment

Product Code	Product Description	Pack Size	Quote
MWC-10-0000006 09	Topic 20 Floor Mounted Bedpan Washer	ea	
MWC-10-0001399 33-A	Topic 20 A Floor Mounted Bedpan Washer with Slop Sink	ea	

Relevant Product

- [Asepti Synergy WD, 5 L](#)
- [Asepti Autozyme Pr, 5 L](#)

- [Asepti Synergy 5, 5L](#)
- [Asepti Wash, 5L](#)
- [Asepti Sonic, 5L](#)
- [Asepti RHW, 5L](#)
- [Asepti Rinse Clean, 5L](#)

Applications

- Infection Control
- Healthcare Facilities
- Residential Aged Care

Resources

- Video
- Brochure