

MI Transporter

Ecolab's Medical Instrument (MI) trolley is the perfect solution for busy hospitals and day surgeries.





ARTG 474328



Features

- Ecolab's medical instrument transportation trolley is the ideal solution for safely transporting your medical instruments or endoscopy equipment throughout your facility.
- MI Transporter feature large drawers to easily accommodate your choice of instruments; 510 x 510 x 100mm
- Made from tough material with a durable, ergonomic design, MI Transporters are built to last.
- MI Transporters feature individually lockable wheel casters, making transport a breeze and provides safe, secure use at bedside.
- Ecolab MI Transporters are available in three different sizes to suit your departments needs; 5 drawer to fit under bench, 6 drawer for standard bench height and 10 drawer ideal for bulk short term storage or transportation of multiple instruments for busy surgeries.
- MI Transporters are built with smooth surfaces for easy cleaning and to assist with hygienic infection control.
- Sterile, single use, disposable liners are available to protect the tray with disposable tray covers to enclose medical instruments to protect and minimise risk of cross contamination.
- Ecolab's MI Transporter tray covers are double sided with easy identification of clean or dirty contents; Green for clean, Red for dirty.
- MI Transporter trays are easy to clean; simply wipe over using an environmental surface cleaner/disinfectant or trays may be reprocessed in a cart washer at low temperature (no higher than 60°C). Refer to product IFU for further instructions and direction.
- Ecolab's MI Transport Trolleys are proudly manufactured in Australia.
- A convenient wall mountable dispenser is available for tray liners and covers for quick, easy and hygienic access.

Product Code	Product Description	Pack Size	Quote
MI-20095, 30000262	Ecolab MI Transport Trolley, 5 Draw	ea	
MI-20096, 30000263	Ecolab MI Transport Trolley, 6 Draw	ea	

Product Code	Product Description	Pack Size	Quote
MI-20010, 30000260	Ecolab MI Transport Trolley, 10 Draw	ea	
MI-20099, 10117894	Dispenser for draw liners/covers	ea	
MI-20091, 30000261	Liners to suit Ecolab MI Transportation Trolley	50 Pk	
MI-C10081, 30000264	Liners to suit Cleanascope Transportation Trolley	50 Pk	

Relevant Product

- [Soluscope Sprint](#)
- [Soluscope Sprint PT](#)
- [Soluscope DSC8000](#)
- [Asepti Active Wipes](#)
- [Asepti Active Liquid](#)
- [Asepti Dry Wipes](#)
- [Belimed WD290IQ](#)

Applications

- Endoscopy
- Healthcare Facilities
- Central Sterile Department

Resources

- Video
- Brochure