


Smartline SmartSink






















Electric Height Adjustable Washing Stations improve productivity










Features

- Smartline SmartSink offers your reprocessing department an electric height adjustable washing station which improves productivity and reduces the risk of staff injury.
- Smartline SmartSinks are manufactured in Australia and constructed out of 304 grade stainless steel. The Smartline range of SmartSinks offer a complete, customisable washing station solution for your departments requirements.
- Smartline SmartSinks feature a 300mm high tapered splashback and tapered benchtop with wetdeck perimeter, offering smooth surfaces for easy cleaning for infection prevention.
- SmartSinks include easy-to-use electric height adjustment with anti-collision technology from 850mm to 1080mm off the ground. Controls are positioned under the sink top offering the user full control and easy access.
- Smartline SmartSinks offer so many customisable options; multiple bench surface, two different or customisable bowl depth, rectangular, circular, single and/or double sink options, sink inserts, shelves, trays and racks and tapware to name a few.
- All Smartline SmartSinks are manufactured with adjustable bolt down feet.

Product Code	Product Description	Pack Size	Quote
SMT-11.SS-135-SC, 10118351	Smartline 1350mm SmartSink Single Bowl - Circular	ea	
SMT-11.SS-175-SC, 10118355	Smartline 1750mm SmartSink Single Bowl - Circular	ea	
sSMT-11.SS-195-SC, 10118359	Smartline 1950mm SmartSink Single Bowl - Circular	ea	
SMT-11.SS-220-SC, 10118364	Smartline 2200mm SmartSink Single Bowl - Circular	ea	

Product Code	Product Description	Pack Size	Quote
SMT-11.SS-175-DC, 10118353	Smartline 1750mm SmartSink Double Bowl - Circular	ea	
SMT-11.SS-195-DC, 10118357	Smartline 1950mm SmartSink Double Bowl - Circular	ea	
SMT-11.SS-220-DC, 10118361	Smartline 2200mm SmartSink Double Bowl - Circular	ea	
SMT-11.SS-260-DC, 10118366	Smartline 2600mm SmartSink Double Bowl - Circular	ea	
SMT-11.SS-320-DC, 10118370	Smartline 3200mm SmartSink Double Bowl - Circular	ea	
SMT-11.SS-360-DC, 10118373	Smartline 3600mm SmartSink Double Bowl - Circular	ea	
SMT-11.SS- UT13560, 10118381	Smartline Under Sink Tray - suit 1350x600 Sink Units	ea	
SMT-11.SS- UT13570, 10118382	Smartline Under Sink Tray - suit 1350x700 Sink Units	ea	
SMT-11.SS- UT17560, 10118383	Smartline Under Sink Tray - suit 1750x600 Sink Units	ea	
SMT-11.SS- UT17570, 10118384	Smartline Under Sink Tray - suit 1750x700 Sink Units	ea	
SMT-11.SS- UT19570, 10118386	Smartline Under Sink Tray - suit 1950x700 Sink Units	ea	
SMT-11.SS- UT22060, 10118387	Smartline Under Sink Tray - suit 2200x600 Sink Units	ea	
SMT-11.SS- UT22070, 10118388	Smartline Under Sink Tray - suit 2200x700 Sink Units	ea	
SMT-11.SS- UT26060, 10118389	Smartline Under Sink Tray - suit 2600x600 Sink Units	ea	
SMT-11.SS- UT26070, 10118390	Smartline Under Sink Tray - suit 2600x700 Sink Units	ea	
SMT-11.SBS2050, 10118346	Smartline Smartsink Splash Back Shelf - 20x50cm	ea	
SMT-11.SBS2080, 10118347	Smartline Smartsink Splash Back Shelf - 20x80cm	ea	
SMT-11.SBT2050, 10118348	Smartline SmartSink Splashback Tray 20x50cm	ea	
SMT-11.SBT2080, 10118349	Smartline SmartSink Splashback Tray 20x80cm	ea	
SMT-11.SS- CHEMSHLF, 10118377	Smartline SmartSink Chemical Bottle Shelf	ea	
SMT-11.SS-SBCM	Smartline SmartSink Clamp Mount for Accessories	ea	

Product Code	Product Description	Pack Size	Quote
SMT-51.TRYRCK-1X5, 10118577	Smartline SmartSink Stainless Steel Tray Rack for Autoplas Trays	ea	
SMT-51.TRDIN-1X5, 10118576	Smartline SmartSink Stainless Steel Tray Rack for DIN Trays	ea	
SMT-11.SS-GF152510-SS-SH, 10118576	Smartline SmartSink Pre-Rinse Unit, 10 inch spout	ea	
SMT-11.TAP-CL10010, 10118393	Smartline SmartSink Lever Mixing Set	ea	
SMT-11.TAP-CL15035, 10118394	Smartline SmartSink Lever Sink Set	ea	
SMT-11.TAP-JETF3000W, 10118397	Smartline SmartSink Pre-Rinse Unit	ea	
SMT-47.DINA-JET, 10118512	Smartline SmartSink Cleaning Jet Pistol with Standard Set	ea	
SMT-47.HELA-JET, 10118513	Smartline SmartSink Cleaning Jet Pistol with Standard Set, Stainless Steel	ea	

Relevant Product

- [Asepti Epizyme Ultra](#)
- [Asepti Multizyme](#)
- [Asepti Neutrazyme](#)
- [Asepti Rapidzyme PR](#)
- [Asepti Release Plus](#)
- [Asepti Synergy 5](#)
- [Asepti Synergy WD](#)
- [Asepti Wash](#)
- [Optipro Multi Enzyme](#)
- [Water Testing, Table 7.3](#)
- [EndoSelect Cleaning Brushes](#)
- [EndoSelect Cleaning Sponges](#)
- [Transportation Trolleys](#)
- [Central Service Program](#)

Applications

- Endoscopy
- Central Sterile Department
- Healthcare Facilities

Resources

- video
- Brochure